

# 1- door tall safety cabinet Type 30 for gas cylinders

## COMPLIANCE

EN 14470-2, EN 1363-1

## ADVANTAGES

- 30-minute fire resistant and made in France
- EN 14470-2 and EN 1363-1 compliant

## BUT ALSO...

- Tested and approved Type 30 according to EN 14470-2 and EN 1363-1
- Normalized labelling provided with the equipment
- Pedestal area under-driveable by a forklift by removing panel
- 30-minute fire resistance
- Double wall construction
- Holes for the connection to the ventilation
- Thermo-expanding seals of the doors
- Adjustable leveling feet for enhance stability
- Locking system



## TECHNICAL FEATURES

Model	<b>7634BG</b>	Door closure	<b>manual closing</b>
Compliance	<b>EN 14470-2, EN 1363-1</b>	Shipping dimensions HxWxD (mm)	<b>2250 x 1200 x 800</b>
Cabinet materials	<b>Built in double wall melamine (thickness: 76mm) :</b>	Ventilation	<b>outlet located at the bottom of the cabinet, inlet located at the top</b>
Doors color	<b>Grey structure (RAL 7035), yellow edges (RAL 1023)</b>	Openings diameter (mm)	<b>100</b>
Number of doors	<b>1</b>	Connecting sleeve	<b>1</b>
Doors type	<b>swing doors</b>	HS Code	<b>246.00</b>
Hinges	<b>right</b>		<b>94032080</b>

## ADDITIONAL ACCESSORIES RECOMMENDED - FILTRATION OR VENTILATION BOXES

- Horizontal rail for mounting 2 cylinders B50 (Réf : RH4BG)
- Galvanized steel ramp for 7634BG or 7635BG (Réf : RABG)

## OPTIONS

- Shelf for 1 cylinder type B2, B5, B10, B11, B15 or B20 (Réf : EG1)

Model	External dimensions HxWxD (mm)	Internal dimensions HxWxD (mm)	Storage capacity (L)	Containment sump	Shelves	Load capacity	Weight (kg)
7634BG	2055 x 660 x 665	1753 x 508 x 458	2		0	0	246.00