



1-door large sheet steel safety cabinet

COMPLIANCE

ADVANTAGES

- Compliant to the Public Health Code
- INRS recommendations compliant
- Storage of flammables, corrosives and toxics

BUT ALSO...

- Compliant to the Public Health Code articles R. 5132-66 and R. 5132-68
- INRS recommendations compliant
- Holes for the connection to the ventilation
- Allows the storage of your chemicals safely
- Normalized labelling provided with the equipment
- Affordable equipment to store your chemicals



TECHNICAL FEATURES

Model	AL155	Shipping dimensions HxWxD (mm)	2000 x 1200 x 800
Cabinet materials	8/10th epoxy painted steel structure	Ventilation	Both inlet and outlet
Doors color	grey RAL 7035	Openings diameter (mm)	100
Number of doors	1	Number of compartments	1
Doors type	swing doors	Containment sump	1
Hinges	right		41.00
Door closure	manual closing	HS Code	94032080

ADDITIONAL ACCESSORIES RECOMMENDED – FILTRATION OR VENTILATION BOXES

- Filtration box with activated carbon filter for organic vapors (Réf : CDF-A)
- Filtration box with activated carbon filter for organic and corrosive vapors (Réf : CDF-ACORG)
- Filtration box built in PVC with activated carbon filter for organic and corrosive vapors (Réf : CDF-P)
- Ventilation box for a connection to the outside (Réf : CDV-A)
- Recirculating air box Labopur - NF X 15-211 (Réf : H50C)
- Connecting kit (Réf : KRC)

OPTIONS

- Perforated cover plate for 1-door cabinet (Réf : CAL155)
- Containment sump for 1-door cabinet (Réf : ESAL155)
- Containment sump in PVC for 1-door cabinet (Réf : SEAL155)

Model	External dimensions HxWxD (mm)	Internal dimensions HxWxD (mm)	Storage capacity (L)	Containment sump	Shelves	Load capacity	Weight (kg)
AL155	1800 x 500 x 420	1730 x 480 x 390	130		3	30	41.00