

APC62+



INRS



2- door tall melamine safety cabinet for acids and bases

COMPLIANCE

ADVANTAGES

- Made in France
- INRS recommendations compliant
- Storage of corrosives

BUT ALSO...

- Adjustable leveling feet for enhance stability
- Adjustable leveling feet for enhance stability
- Cabinet with compartments to separate acids and bases
- PP containment sumps on sliding rails
- INRS recommendations compliant
- Built in melamine with a great resistance to corrosive vapors
- Holes for the connection to the ventilation
- Normalized labelling provided with the equipment
- Especially designed to store acids and bases



TECHNICAL FEATURES

Model	APC62+	Shipping dimensions HxWxD (mm)	870 x 600 x 800
Cabinet materials	Melamine wood structure	Ventilation	Both inlet and outlet
Doors color	white RAL 9010	Openings diameter (mm)	100
Number of doors	2	Number of compartments	2
Doors type	sliding doors		60.00
Hinges	right & left	HS Code	94036000
Door closure	manual closing		

ADDITIONAL ACCESSORIES RECOMMENDED – FILTRATION OR VENTILATION BOXES

- Filtration box with activated carbon filter for organic and corrosive vapors (Réf : CDF-ACORG)
- Filtration box built in PVC with activated carbon filter for organic and corrosive vapors (Réf : CDF-P)
- Connecting kit (Réf : KRC)
- Vertical recycling air box (delivered without filter) (Réf : H40+)

OPTIONS

- Containment sump built in polypropylene (Réf : E8A2)

Model	External dimensions HxWxD (mm)	Internal dimensions HxWxD (mm)	Storage capacity (L)	Containment sump	Shelves	Load capacity	Weight (kg)
APC62+	630 x 1135 x 578		72		0	20	60.00