



## 2- door underbench safety cabinet for acids and flammables

### COMPLIANCE

### ADVANTAGES

- INRS recommendations compliant
- Designed to be placed under-bench
- Storage of flammables, corrosives and toxics

### BUT ALSO...

- Especially designed to settle under the bench
- Equipment allowing space saving
- INRS recommendations compliant
- Holes for the connection to the ventilation
- Normalized labelling provided with the equipment
- Doors with key locking system
- Possibility to add swivel castors, for an easy movement of the cabinet



### TECHNICAL FEATURES

|                   |                                              |                                |                              |
|-------------------|----------------------------------------------|--------------------------------|------------------------------|
| Model             | <b>AS60</b>                                  | Shipping dimensions HxWxD (mm) | <b>870 x 1200 x 800</b>      |
| Cabinet materials | <b>15/10th epoxy painted steel structure</b> | Ventilation                    | <b>Both inlet and outlet</b> |
| Doors color       | <b>blue RAL 5015 &amp; yellow RAL 1024</b>   | Openings diameter (mm)         | <b>100</b>                   |
| Number of doors   | <b>2</b>                                     | Number of compartments         | <b>2</b>                     |
| Doors type        | <b>swing doors</b>                           | Containment sump               | <b>60.00</b>                 |
| Hinges            | <b>right &amp; left</b>                      | HS Code                        | <b>94032080</b>              |
| Door closure      | <b>self-closing</b>                          |                                |                              |

### ADDITIONAL ACCESSORIES RECOMMENDED - FILTRATION OR VENTILATION BOXES

- Filtration box with activated carbon filter for organic and corrosive vapors (Réf : CDF-ACORG)
- Filtration box built in PVC with activated carbon filter for organic and corrosive vapors (Réf : CDF-P)
- Ventilation box for a connection to the outside (Réf : CDV-A)
- Connecting kit (Réf : KRC)
- Vertical recycling air box (delivered without filter) (Réf : H40+)

### OPTIONS

- 4 swivel castors 2 with locking (Réf : 4ROU)

| Model | External dimensions HxWxD (mm) | Internal dimensions HxWxD (mm) | Storage capacity (L) | Containment sump | Shelves | Load capacity | Weight (kg) |
|-------|--------------------------------|--------------------------------|----------------------|------------------|---------|---------------|-------------|
| AS60  | 670 x 1100 x 500               | (2x) 650 x 530 x 470           | 60                   |                  | 2       | 50            | 60.00       |