



1- red door working cover safety cabinet

COMPLIANCE

ADVANTAGES

- Compliant to the Public Health Code
- INRS recommendations compliant
- Storage of flammables, corrosives and toxics

BUT ALSO...

- Compliant to the Public Health Code articles R. 5132-66 and R. 5132-68
- INRS recommendations compliant
- Allows the storage of your chemicals safely
- Normalized labelling provided with the equipment
- Affordable equipment to store your chemicals



TECHNICAL FEATURES

Model	AZ150R	Shipping dimensions HxWxD (mm)	2150 x 1200 x 800
Cabinet materials	8/10th epoxy painted steel structure	Ventilation	Both inlet and outlet
Doors color	red RAL 3000	Number of compartments	1
Number of doors	1	Containment sump	1
Doors type	swing doors	HS Code	55.00
Hinges	right		94032080
Door closure	manual closing		

ADDITIONAL ACCESSORIES RECOMMENDED – FILTRATION OR VENTILATION BOXES

- Opening vent for range 14.E (Réf : AZVENT)
- Filtration box with activated carbon filter for organic vapors (Réf : CDF-A)
- Filtration box with activated carbon filter for organic and corrosive vapors (Réf : CDF-ACORG)
- Filtration box built in PVC with activated carbon filter for organic and corrosive vapors (Réf : CDF-P)
- Ventilation box for a connection to the outside (Réf : CDV-A)
- Recirculating air box Labopur - NF X 15-211 (Réf : H50C)
- Connecting kit (Réf : KRC)

OPTIONS

- Perforated cover plate for cabinets AZ50 - AZ150 - AC150 (Réf : CE150)
- Containment sump for AZ50 - AZ150 - AC150 (Réf : ESE150)

Model	External dimensions HxWxD (mm)	Internal dimensions HxWxD (mm)	Storage capacity (L)	Containment sump	Shelves	Load capacity	Weight (kg)
AZ150R	1950 x 500 x 500	1840 x 498 x 453	150		3	80	55.00