



1- door underbench safety cabinet for flammables

COMPLIANCE

ADVANTAGES

- INRS recommendations compliant
- Designed to be placed under-bench
- Storage of flammables, corrosives and toxics

BUT ALSO...

- Equipment allowing space saving
- Especially designed to settle under the bench
- Equipment allowing space saving
- Normalized labelling provided with the equipment
- INRS recommendations compliant
- Doors with key locking system
- Possibility to add swivel castors, for an easy movement of the cabinet
- Holes for the connection to the ventilation



TECHNICAL FEATURES

Model	SV30	Shipping dimensions HxWxD (mm)	870 x 1200 x 800
Cabinet materials	15/10th epoxy painted steel structure	Ventilation	Both inlet and outlet
Doors color	yellow RAL 1023	Openings diameter (mm)	100
Number of doors	1	Number of compartments	1
Doors type	swing doors	Containment sump	1
Hinges	left		30.00
Door closure	self-closing	HS Code	94032080

ADDITIONAL ACCESSORIES RECOMMENDED – FILTRATION OR VENTILATION BOXES

- Filtration box with activated carbon filter for organic vapors (Réf : CDF-A)
- Ventilation box for a connection to the outside (Réf : CDV-A)
- Connecting kit (Réf : KRC)
- Vertical recycling air box (delivered without filter) (Réf : H40+)

OPTIONS

- 4 swivel castors 2 with locking (Réf : 4ROU)

Model	External dimensions HxWxD (mm)	Internal dimensions HxWxD (mm)	Storage capacity (L)	Containment sump	Shelves	Load capacity	Weight (kg)
SV30	670 x 550 x 500	650 x 530 x 470	30		1	50	30.00